

Phospho-Tau (Ser404) Rabbit mAb

ABZ4524


Product Overview

Name	Phospho-Tau (Ser404) Rabbit mAb
Accession(Primary)	P10636
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 79 kDa; Observed MW: 50-80 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:20-1:50; IP-1:20

Target information (P10636)

Gene Name	MAPT
About Gene Symbol	

Human Fetal brain



Western blot analysis of Tau (Phospho-Ser404) in Human Fetal brain lysates using Phospho-Tau (Ser404) antibody.