



Phospho-ATF2 (Thr71) Rabbit mAb

ABZ4526


Product Overview

Name	Phospho-ATF2 (Thr71) Rabbit mAb
Accession(Primary)	P15336
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 55 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P15336)

Gene Name ATF2**About Gene Symbol**

Immunocytochemistry analysis of Phospho-ATF2 (Thr71) (green) in U87-MG using Phospho-ATF2 (Thr71) antibody, and DAPI(blue).



Western blot analysis of Phospho-ATF2 (Thr71) in 3T3, HeLa lysates using Phospho-ATF2 (Thr71) antibody.