

Phospho-Cannabinoid Receptor I (Ser316) Rabbit mAb


ABZ4528

Product Overview

Name	Phospho-Cannabinoid Receptor I (Ser316) Rabbit mAb
Accession(Primary)	P21554
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P21554)

Gene Name	CNR1
About Gene Symbol	



Western blot analysis of Phospho-Cannabinoid Receptor I (Ser316) in HeLa lysates using Phospho-Cannabinoid Receptor I (Ser316) antibody.