

Phospho-CBL (Ser669) Rabbit mAb

ABZ4529


Product Overview

Name	Phospho-CBL (Ser669) Rabbit mAb
Accession(Primary)	P22681
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 100 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P22681)

Gene Name CBL

About Gene Symbol



Western blot analysis of Phospho-CBL (Ser669) in K562 lysates using Phospho-CBL (Ser669) antibody.

