


Phospho-CD18 (Ser745) Rabbit mAb

ABZ4530

Product Overview

Name	Phospho-CD18 (Ser745) Rabbit mAb
Accession(Primary)	P05107
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:50-1:100

Target information (P05107)

Gene Name ITGB2**About Gene Symbol**

Western blot analysis of CD18 (Phospho-Ser745) in K562, rat Brain, C6, 3T3 lysates using Phospho-CD18 (Ser745) antibody.

