

Phospho-CD3 zeta (Tyr142) Rabbit mAb


ABZ4531

Product Overview

Name	Phospho-CD3 zeta (Tyr142) Rabbit mAb
Accession(Primary)	P20963
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 19 kDa; Observed MW: 19 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (P20963)

Gene Name	CD247
About Gene Symbol	



Western blot analysis of Phospho-CD3 zeta (Tyr142) in Jurkat lysates using Phospho-CD3 zeta (Tyr142) antibody.