

## Phospho-CDK1/2/3 (Thr14) Rabbit mAb


ABZ4532

### Product Overview

|                             |                                                                         |
|-----------------------------|-------------------------------------------------------------------------|
| <b>Name</b>                 | Phospho-CDK1/2/3 (Thr14) Rabbit mAb                                     |
| <b>Accession(Primary)</b>   | P06493                                                                  |
| <b>Host Species</b>         | Rabbit                                                                  |
| <b>Clonality</b>            | Monoclonal                                                              |
| <b>Isotype</b>              | IgG                                                                     |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                       |
| <b>Calculated MW</b>        | Calculated MW: 34 kDa; Observed MW: 34 kDa                              |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:50-1:200;<br>ICC/IF-1:100-1:500; IP-1:20-1:50 |

### Target information (P06493)

|                          |                |
|--------------------------|----------------|
| <b>Gene Name</b>         | CDK1/CDK2/CDK3 |
| <b>About Gene Symbol</b> |                |



Western blot analysis of Phospho-  
CDK1/CDK2/CDK3 (Thr14) in HeLa, CHO-K1 lysates  
using Phospho-CDK1/2/3 (Thr14) antibody.