

Phospho-Cytokeratin 8 (Ser23) Rabbit mAb


ABZ4534

Product Overview


Name	Phospho-Cytokeratin 8 (Ser23) Rabbit mAb
Accession(Primary)	P05787
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:50-1:100

Target information (P05787)


Gene Name	KRT8
About Gene Symbol	



Western blot analysis of Phospho-Cytokeratin 8 (Ser23) in K562, HeLa lysates using Phospho-Cytokeratin 8 (Ser23) antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Cytokeratin 8 (Phospho-Ser23) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Cytokeratin 8
(Phospho-Ser23)(green) in HeLa using Cytokeratin 8
(Phospho-Ser23) antibody, and DAPI(blue)