

Phospho-Dishevelled 2 (Ser143) Rabbit mAb

ABZ4535

Product Overview

Name	Phospho-Dishevelled 2 (Ser143) Rabbit mAb
Accession(Primary)	O14641
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 79 kDa; Observed MW: 90-95 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O14641)

Gene Name DVL2

About Gene Symbol



