

Phospho-IGF1 Receptor (Tyr1166) Rabbit mAb

ABZ4538


Product Overview

Name	Phospho-IGF1 Receptor (Tyr1166) Rabbit mAb
Accession(Primary)	P08069
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 155 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (P08069)

Gene Name IGF1R

About Gene Symbol



Western blot analysis of Phospho-IGF1R (Tyr1166) in C6, 3T3, HeLa lysates using Phospho-IGF1 Receptor (Tyr1166) antibody.

