

Phospho-LAT (Tyr200) Rabbit mAb


ABZ4539

Product Overview

Name	Phospho-LAT (Tyr200) Rabbit mAb
Accession(Primary)	O43561
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 28 kDa; Observed MW: 36-38 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:200-1:2000

Target information (O43561)

Gene Name	LAT
About Gene Symbol	LAT (Linker For Activation Of T Cells) is a signaling adaptor protein involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-LAT (Tyr200) in Jurkat lysates using Phospho-LAT (Tyr200) antibody.