

Phospho-Nucleolin (Thr84) Rabbit mAb

ABZ4540


Product Overview

Name	Phospho-Nucleolin (Thr84) Rabbit mAb
Accession(Primary)	P19338
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 77 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000

Target information (P19338)

Gene Name NCL

About Gene Symbol



Western blot analysis of Phospho-Nucleolin (Thr84) in U251MG, HepG2, Hela lysates using Phospho-Nucleolin (Thr84) antibody.

