

Phospho-Parvin alpha (Ser8) Rabbit mAb

ABZ4541


Product Overview

Name	Phospho-Parvin alpha (Ser8) Rabbit mAb
Accession(Primary)	Q9NVD7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200

Target information (Q9NVD7)

Gene Name PARVA

About Gene Symbol



Western blot analysis of Parvin alpha (Phospho-Ser8) in U2OS, mouse brain, rat brain, PC-12, A549 lysates using Phospho-Parvin alpha (Ser8) antibody.

