

Phospho-PRC1 (Thr481) Rabbit mAb


ABZ4543

Product Overview


Name	Phospho-PRC1 (Thr481) Rabbit mAb
Accession(Primary)	O43663
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 72 kDa; Observed MW: 72 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; ICC/IF-1:50-1:100

Target information (O43663)


Gene Name	PRC1
About Gene Symbol	PRC1 (Protein Regulator Of Cytokinesis 1) is a cytokine and immune signaling protein involved in JAK-STAT/cytokine signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-PRC1 (Thr481) in HeLa, Jurkat lysates using Phospho-PRC1 (Thr481) antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using PRC1 (Phospho-Thr481) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Phospho-PRC1
(Thr481) (green) in HeLa using Phospho-PRC1
(Thr481) antibody, and DAPI (blue)