

Phospho-Progesterone Receptor (Ser190) Rabbit mAb


ABZ4544

Product Overview

Name	Phospho-Progesterone Receptor (Ser190) Rabbit mAb
Accession(Primary)	P06401
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 99 kDa; Observed MW: 90,118 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (P06401)

Gene Name	PGR
About Gene Symbol	



Western blot analysis of Progesterone Receptor (Phospho-Ser190) in T47D lysates using Phospho- Progesterone Receptor (Ser190) antibody.