

Phospho-Rb2 p130 (Ser952) Rabbit mAb


ABZ4547

Product Overview


Name	Phospho-Rb2 p130 (Ser952) Rabbit mAb
Accession(Primary)	Q08999
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 128 kDa; Observed MW: 128 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; ICC/IF-1:50-1:100

Target information (Q08999)


Gene Name	RBL2
About Gene Symbol	



Western blot analysis of Rb2 p130 (Phospho-Ser952) in K562, HeLa lysates using Phospho-Rb2 p130 (Ser952) antibody.



Immunocytochemistry analysis of Rb2 p130 (Phospho-Ser952)(green) in HeLa using Rb2 p130 (Phospho-Ser952) antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Rb2 p130 (Phospho-Ser952) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.