

Phospho-Tau (Ser202) Rabbit mAb

ABZ4548


Product Overview

Name	Phospho-Tau (Ser202) Rabbit mAb
Accession(Primary)	P10636
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 79 kDa; Observed MW: 50-80 kDa
Recommended Dilution	WB-1:1000-1:5000


Target information (P10636)

Gene Name MAPT

About Gene Symbol



Western blot analysis of Phospho-Tau (Ser202) in mouse brain, HeLa, HeLa+CalyculinA lysates using Phospho-Tau (Ser202) antibody.



Western blot analysis of Phospho-Tau (Ser202) in rat Brain lysates using Phospho-Tau (Ser202) antibody.