

Phospho-Tau (Thr50) Rabbit mAb


ABZ4550

Product Overview


Name	Phospho-Tau (Thr50) Rabbit mAb
Accession(Primary)	P10636
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 79 kDa; Observed MW: 50-80 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20; IHC-1:50-1:100

Target information (P10636)

Gene Name	MAPT
About Gene Symbol	



Immunocytochemistry analysis of Phospho-Tau (Thr50) (green) in HCT116 using Phospho-Tau (Thr50) antibody, and DAPI (blue).



Western blot analysis of Tau (Phospho-Thr50) in Human Hippocampus lysates using Tau (Phospho-Thr50) antibody.