

Phospho-beta Arrestin 1 (Ser412) Rabbit mAb


ABZ4554

Product Overview


Name	Phospho-beta Arrestin 1 (Ser412) Rabbit mAb
Accession(Primary)	P49407
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 47 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:50-1:100

Target information (P49407)

Gene Name	ARRB1
About Gene Symbol	ARRB1 (Arrestin Beta 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-beta Arrestin 1 (Ser412) in A549, MCF-7, U251, HL-60 lysates using Phospho-beta Arrestin 1 (Ser412) antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using beta Arrestin 1 (Phospho-Ser412) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.