

Phospho-EEF2/Elongation factor 2 (Thr56/Thr58) Rabbit mAb


ABZ4560

Product Overview

Name	Phospho-EEF2/Elongation factor 2 (Thr56/Thr58) Rabbit mAb
Accession(Primary)	P13639
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 95 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P13639)

Gene Name	EEF2
About Gene Symbol	EEF2 (Eukaryotic Translation Elongation Factor 2) is a metabolic enzyme involved in translational regulation. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of EEF2/Elongation factor 2 (Phospho-Thr56/Thr58) in HeLa, HL-60 lysates using Phospho-Elongation Factor 2 (Thr56/Thr58) antibody.