

Phospho-GATA3 (Ser308) Rabbit mAb


ABZ4564

Product Overview


Name	Phospho-GATA3 (Ser308) Rabbit mAb
Accession(Primary)	P23771
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20

Target information (P23771)

Gene Name	GATA3
About Gene Symbol	GATA3 (GATA Binding Protein 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-GATA3 (Ser308) in A549, MCF-7, U251 lysates using Phospho-GATA3 (Ser308) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using GATA3 (Phospho-Ser308) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.