

Phospho-Hormone Sensitive Lipase/HSL (Ser855) Rabbit


ABZ4565

Product Overview

Name	Phospho-Hormone Sensitive Lipase/HSL (Ser855) Rabbit mAb
Accession(Primary)	Q05469
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 117 kDa; Observed MW: 83 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q05469)

Gene Name	LIPE
About Gene Symbol	



Western blot analysis of Hormone sensitive lipase/HSL (Phospho-Ser855) in A549, HL-60, U2OS lysates using Phospho-Hormone Sensitive Lipase (Ser855) antibody.