

Phospho-c-Jun/JunD (Ser73/Ser100) Rabbit mAb


ABZ4566

Product Overview

Name	Phospho-c-Jun/JunD (Ser73/Ser100) Rabbit mAb
Accession(Primary)	P05412
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 38,42,48 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; IP-1:50-1:100

Target information (P05412)

Gene Name	JUN/JUND
About Gene Symbol	JUN (Jun Proto-Oncogene, AP-1 Transcription Factor Subunit) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Jun (Phospho-Ser73)/JunD (Phospho-Ser100) in A549, HL-60, U251, U87-MG lysates using Phospho-Jun/JunD (Ser73/Ser100) antibody.