

Phospho-c-Met (Tyr1349) Rabbit mAb


ABZ4567

Product Overview

Name	Phospho-c-Met (Tyr1349) Rabbit mAb
Accession(Primary)	P08581
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 156 kDa; Observed MW: 170,140 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; IP-1:20-1:200

Target information (P08581)

Gene Name	MET
About Gene Symbol	



Western blot analysis of Phospho-Met (c-Met) (Tyr1349) in HeLa, A549 lysates using Phospho-c-Met (Tyr1349) antibody.