

Phospho-Nucleophosmin (Thr199) Rabbit mAb


ABZ4569

Product Overview


Name	Phospho-Nucleophosmin (Thr199) Rabbit mAb
Accession(Primary)	P06748
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 33 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20; IF-1:50-1:200

Target information (P06748)

Gene Name	NPM1
About Gene Symbol	



Western blot analysis of Nucleophosmin (Phospho-Thr199) in C6, Ramos lysates using Phospho-Nucleophosmin (Thr199) antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Nucleophosmin (Phospho-Thr199) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.