

Phospho-PTP1B (Ser378) Rabbit mAb

ABZ4572


Product Overview


Name	Phospho-PTP1B (Ser378) Rabbit mAb
Accession(Primary)	P18031
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (P18031)

Gene Name PTPN1

About Gene Symbol





Immunohistochemistry analysis of paraffin-embedded Human tonsil using PTP1B (Phospho-Ser378) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.