

Phospho-Rad17 (Ser656) Rabbit mAb


ABZ4574

Product Overview


Name	Phospho-Rad17 (Ser656) Rabbit mAb
Accession(Primary)	O75943
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 77 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; IHC-P-1:100-1:200

Target information (O75943)


Gene Name	RAD17
About Gene Symbol	RAD17 (RAD17 Checkpoint Clamp Loader Component) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Phospho-Rad17 (Ser656) (green) in HeLa using Phospho-Rad17 (Ser656) antibody, and DAPI (blue).



Western blot analysis of Rad17 (Phospho-Ser656) in HeLa, Ramos lysates using Rad17 (Phospho-Ser656) antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using Rad17 (Phospho-Ser656) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.