

Phospho-Src (Tyr529) Rabbit mAb


ABZ4579

Product Overview


Name	Phospho-Src (Tyr529) Rabbit mAb
Accession(Primary)	P12931
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:500

Target information (P12931)

Gene Name	SRC
About Gene Symbol	SRC (SRC Proto-Oncogene, Non-Receptor Tyrosine Kinase) is a receptor tyrosine kinase involved in receptor-mediated signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-Src (Tyr529) in rat lung lysates using Phospho-Src (Tyr530) antibody.



Western blot analysis of Src (Phospho-Tyr529) in CHO-K1 lysates using Src (Phospho-Tyr529) antibody.