

Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) Rabbit mAb


ABZ4580

Product Overview

Name	Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) Rabbit mAb
Accession(Primary)	O96013
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 64,80,75 kDa; Observed MW: 64,80,75 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O96013)

Gene Name	PAK4/PAK5/PAK6
About Gene Symbol	



Western blot analysis of Phospho-PAK4 + PAK5 + PAK6 (Ser474 + Ser560 + Ser602) in HeLa, A549, HL-60, U2OS, U87-MG lysates using Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) antibody.