

Phospho-CREB (Ser133) Rabbit pAb


ABZ4581

Product Overview

Name	Phospho-CREB (Ser133) Rabbit pAb
Accession(Primary)	P16220
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Pig
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; IP-2-5 ug/mg lysate; ELISA-1:10000

Target information (P16220)

Gene Name	CREB1
About Gene Symbol	CREB1 (CAMP Responsive Element Binding Protein 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular



Western blot analysis of Phospho-CREB (Ser133)
Phospho-CREB-1 (S133) in various lysates using Phospho-CREB (Ser133)
antibody.