

Phospho-Myosin Phosphatase (Thr853) Rabbit pAb

ABZ4588

Product Overview

Name	Phospho-Myosin Phosphatase (Thr853) Rabbit pAb
Accession(Primary)	O14974
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 115 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (O14974)

Gene Name	PPP1R12A
About Gene Symbol	