

Phospho-GARNL1 (Thr736) Rabbit pAb

ABZ4589

Product Overview

Name	Phospho-GARNL1 (Thr736) Rabbit pAb
Accession(Primary)	Q6GYQ0
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 230 kDa; Observed MW: 220 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:200; IF-1:50-1:200

Target information (Q6GYQ0)

Gene Name	RALGAPA1
About Gene Symbol	