


Phospho-alpha Tubulin (Tyr272) Rabbit mAb

ABZ4591

Product Overview

Name	Phospho-alpha Tubulin (Tyr272) Rabbit mAb
Accession(Primary)	P68363
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 50 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (P68363)

Gene Name TUBA1B**About Gene Symbol**

Western blot analysis on Jurkat lysates treated with pervanadate using Phospho-alpha Tubulin (Tyr272) antibody.

