

Phospho-p53 (Ser392) Rabbit mAb


ABZ4592

Product Overview

Name	Phospho-p53 (Ser392) Rabbit mAb
Accession(Primary)	P04637
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:44 kDa;Observed MW: 53 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:20

Target information (P04637)

Gene Name	TP53
About Gene Symbol	TP53 (Tumor Protein P53) is a tumor suppressor protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-p53 (Ser392) in HEK293, A431 lysates using Phospho-p53 (Ser392) antibody.