

Phospho-Rb (Ser807) Rabbit mAb


ABZ4593

Product Overview

Name	Phospho-Rb (Ser807) Rabbit mAb
Accession(Primary)	P06400
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:106 kDa;Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IF-1:200-1:500

Target information (P06400)

Gene Name	RB1
About Gene Symbol	RB1 (RB Transcriptional Corepressor 1) is a protein phosphatase involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-Rb (Ser807) in K562, C6 lysates using Phospho-Rb (Ser807) antibody.