

Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) Rabbit mAb


ABZ4595

Product Overview

Name	Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) Rabbit mAb
Accession(Primary)	O96013
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 64,80,75 kDa;Observed MW:64kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:50


Target information (O96013)

Gene Name	PAK4/5/6
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of Phospho-PAK4/5/6
(Ser474/Ser560/Ser602) in Hela lysates using
Phospho-PAK4/5/6 (Ser474/Ser560/Ser602)
antibody.



Western blot analysis of Phospho-PAK4 + PAK5 +
PAK6 (S474 + S560 + S602) in U87-MG lysates
using Phospho-PAK4/5/6 (Ser474/Ser560/Ser602)
antibody.