

Phospho-Pyruvate Dehydrogenase E1 (Ser293) Rabbit mAb


ABZ4596

Product Overview

Name	Phospho-Pyruvate Dehydrogenase E1 (Ser293) Rabbit mAb
Accession(Primary)	P08559
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:44kDa;Observed MW:40kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:10-1:100


Target information (P08559)

Gene Name	PDHA1
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of Phospho-PDHA1 (Ser293) in 293T lysates using Phospho-PDHA1 (Ser293) antibody.



Western blot analysis of Phospho-PDHA1 (S293) in 293T lysates using Phospho-PDHA1 (Ser293) antibody.