

Lamin A/C Rabbit mAb

ABZ4597



Product Overview

Name	Lamin A/C Rabbit mAb
Accession(Primary)	P02545
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 74 kDa; Observed MW: 74,63 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:100-1:200; F:50-1:100


Target information (P02545)

Gene Name LMNA


About Gene Symbol




Immunocytochemistry analysis of Lamin A/C (green) in HeLa using Lamin A/C antibody, and DAPI(blue).



Western blot analysis of Lamin A/C in 3T3, C2C12 lysates using Lamin A/C antibody.



Western blot analysis of Lamin A/C in U251, Lncap, A549 lysates using Lamin A/C antibody



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Lamin A/C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.