

Lamin A/C (5D12) Mouse mAb

ABZ4598


Product Overview


Name	Lamin A/C (5D12) Mouse mAb
Accession(Primary)	P02545
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 74 kDa; Observed MW: 63,74 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P02545)

Gene Name LMNA

About Gene Symbol





Immunofluorescence analysis of Lamin A/C (5D12)
in A549 using Lamin A/C antibody.