


## Lamin A/C Mouse mAb

ABZ4599

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Lamin A/C Mouse mAb                        |
| <b>Accession(Primary)</b>   | P02545                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG2a                                      |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 74 kDa; Observed MW: 74 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

### Target information (P02545)

**Gene Name** LMNA**About Gene Symbol**

Western blot analysis of Lamin A in K562, C6, 3T3 and HeLa lysates using Lamin A antibody.

