

# TriMethyl-Histone H3 (Lys14) Rabbit pAb


ABZ4600

## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | TriMethyl-Histone H3 (Lys14) Rabbit pAb    |
| <b>Accession(Primary)</b>   | P68431                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Polyclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Rat                                 |
| <b>Calculated MW</b>        | Calculated MW: 16 kDa; Observed MW: 16 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:2000                           |

## Target information (P68431)

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | H3-4                                                                                                                                                                                                                                                            |
| <b>About Gene Symbol</b> | H3C1 (H3 Clustered Histone 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to el |



Western blot analysis of TriMethyl-Histone H3 (Lys14) in various cell lines lysates using TriMethyl-Histone H3K14 antibody.