

MonoMethyl-Histone H3 (Lys4) Rabbit pAb


ABZ4601

Product Overview

Name	MonoMethyl-Histone H3 (Lys4) Rabbit pAb
Accession(Primary)	Q16695
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:500-1:1000

Target information (Q16695)

Gene Name	H3-4
About Gene Symbol	



Western blot analysis of MonoMethyl-Histone H3 (Lys4) in various cell lines lysates using MonoMethyl-Histone H3-K4 antibody.