

TriMethyl-Histone H3 (Lys79) (6A6) Mouse mAb


ABZ4602

Product Overview


Name	TriMethyl-Histone H3 (Lys79) (6A6) Mouse mAb
Accession(Primary)	P68431
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (P68431)


Gene Name	H3C1
About Gene Symbol	H3C1 (H3 Clustered Histone 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to el



Western blot analysis of TriMethyl-Histone H3 (Lys79) in HeLa, rat Testis , Raw264.7 lysates using Histone H3 (tri methyl K79) (3G3) antibody.



Immunoprecipitation analysis of TriMethyl-Histone H3 (Lys79) (6A6) in HeLa lysates using TriMethyl-Histone H3 (Lys79) (6A6) antibody



Western blot analysis of TriMethyl-Histone H3 (Lys79) (6A6) in T47D, COS7 and C6 lysates using TriMethyl-Histone H3 (Lys79) antibody.