

DiMethyl-Histone H3 (Lys27) (5F6) Mouse mAb


ABZ4609

Product Overview


Name	DiMethyl-Histone H3 (Lys27) (5F6) Mouse mAb
Accession(Primary)	P68431
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P68431)


Gene Name	H3C1
About Gene Symbol	H3C1 (H3 Clustered Histone 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to el




Immunoprecipitation analysis of DiMethyl-Histone H3 (Lys27) (5F6) in HeLa lysates using DiMethyl-Histone H3 (Lys27) (5F6) antibody.



Western blot analysis of DiMethyl-Histone H3 (Lys27) (5F6) in C6, 3T3, C2C12 lysates using DiMethyl-Histone H3 (Lys27) (5F6) antibody



Immunocytochemistry analysis of DiMethyl-Histone H3 (Lys27) in HeLa cells using DiMethyl-Histone H3 (Lys27) (5F6) antibody.



Western blot analysis of DiMethyl-Histone H3 (Lys27) (5F6) in HeLa, rat heart, mouse heart, HUVEC, Jurkat lysates using DiMethyl-Histone H3 (Lys27) (5F6) antibody.