

MonoMethyl-Histone H3 (Lys9) (10F5) Mouse mAb


ABZ4611

Product Overview


Name	MonoMethyl-Histone H3 (Lys9) (10F5) Mouse mAb
Accession(Primary)	P68431
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P68431)

Gene Name	H3C1
About Gene Symbol	H3C1 (H3 Clustered Histone 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to el



Western blot analysis of MonoMethyl-Histone H3 in MCF-7, C2C12, C6, 3T3 lysates using MonoMethyl-Histone H3 (Lys9) (10F5) antibody.



Western blot analysis of MonoMethyl-Histone H3 (Lys9) (10F5) in CHO-K1, rat brain, C6, 3T3, HeLa lysates using MonoMethyl-Histone H3 (Lys9) (10F5) antibody.