

MonoMethyl-Histone H3 (Lys18) Rabbit mAb


ABZ4616

Product Overview


Name	MonoMethyl-Histone H3 (Lys18) Rabbit mAb
Accession(Primary)	P68431
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; ICC/IF-1:500-1:5000

Target information (P68431)

Gene Name	H3C1
About Gene Symbol	H3C1 (H3 Clustered Histone 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to el



Western blot analysis of MonoMethyl-Histone H3 (Lys18) in 3T3, rat brain lysates using MonoMethyl-Histone H3 (Lys18) antibody.



Western blot analysis of MonoMethyl-Histone H3 (Lys18) in HeLa, A549, HL-60, U2OS, U87-MG lysates using MonoMethyl-Histone H3 (Lys18) antibody.