

MonoMethyl-Histone H4 (Lys20) Rabbit mAb

ABZ4623


Product Overview

Name	MonoMethyl-Histone H4 (Lys20) Rabbit mAb
Accession(Primary)	P62805
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:11 kDa;Observed MW: 11 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (P62805)

Gene Name H4C1

About Gene Symbol



Western blot analysis of MonoMethyl-Histone H4 (Lys20) in 3T3 lysates using MonoMethyl-Histone H4 (Lys20) antibody.

