

Mono/Di/TriMethyl-Histone H3 (Lys79) Rabbit mAb

ABZ4626

Product Overview


Name	Mono/Di/TriMethyl-Histone H3 (Lys79) Rabbit mAb
Accession(Primary)	P68431
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:15 kDa;Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:500-1:2000

Target information (P68431)

Gene Name

About Gene Symbol

H3C1 (H3 Clustered Histone 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to el



Western blot analysis of Mono/Di/TriMethyl-Histone H3 (Lys79) in HeLa, 3T3 lysates using Mono/Di/TriMethyl-Histone H3 (Lys79) antibody.