


MonoMethyl-Histone H2B (Arg79) Rabbit mAb

ABZ4628

Product Overview

Name	MonoMethyl-Histone H2B (Arg79) Rabbit mAb
Accession(Primary)	P33778
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW:14 kDa;Observed MW: 14 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (P33778)

Gene Name H2BC21**About Gene Symbol**

Western blot analysis of MonoMethyl-Histone H2B (Arg79) in HeLa lysates using MonoMethyl-Histone H2B (Arg79) antibody.

