

## beta I Tubulin Mouse mAb

ABZ4630

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | beta I Tubulin Mouse mAb                   |
| <b>Accession(Primary)</b>   | Q9H4B7                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG1                                       |
| <b>Species Reactivity</b>   | Zebrafish                                  |
| <b>Calculated MW</b>        | Calculated MW: 50 kDa; Observed MW: 50 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

### Target information (Q9H4B7)

**Gene Name** TUBB1**About Gene Symbol**

